

Improved Immobiline™ DryStrip gels and IPG Buffer

Immobiline DryStrip gels (IPG strips) are used for isoelectric focusing (IEF; Fig 1), run as the first dimension of 2-D electrophoresis or as a separate application. IPG strips contain a preformed pH gradient immobilized in precast polyacrylamide gels on a plastic backing. A comprehensive range of overlapping IPG strips covering narrow, medium, and wide pH ranges are available in various strip lengths, and the strips are individually packaged and bar coded for convenient handling and sample tracking. When used for 2-D DIGE applications, the gels give accurate and reliable 2-D results.

Key features:

- Immobiline DryStrip gels have been improved on several levels, through new process parameters, to produce first and second dimension gels of overall higher quality.
- Improved Immobiline DryStrip gels provide increased stability, reproducibility, and protein spot definition.
- Better spot definition improves spot quantitation in data analysis.
- Stable and accurate immobilized pH gradients and precise lengths ensure high reproducibility and reliable gel comparisons.
- Wide selection of IPG strips, pH gradients, lengths, and buffers for greater versatility.
- Immobiline DryStrip gels with matching IPG Buffer gives a high sample loading capacity and improves resolution.

Improved Immobiline DryStrip gels provide increased spot definition

Immobiline DryStrip gels have been improved to increase stability, reproducibility, and protein spot definition. Figure 2A shows a case of skewed bands in the first dimension giving



Fig 1. Using Immobiline DryStrip gels with IPG Buffer maximizes resolution, loading capacity, and reproducibility of first dimension IEF, which greatly improves the quality of your 2-D electrophoresis results. Each IPG strip is individually packaged for extended durability.

rise to elongated spots in the second dimension. In Figure 2B, the improved stability of the pH gradient gives straight bands in the first dimension and better spot definition and spot focusing in the second dimension. With improved Immobiline DryStrip gels, increased spot definition leads to accurate matching and analysis.




Fig 2. Immobiline DryStrip pH 4-7 run according to standard procedures. **A)** Spots are elongated in the second dimension due to skewed bands in the first dimension. **B)** Improved Immobiline DryStrip gels provide a more stable pH gradient and better spot definition in the second dimension.






Fig 3. 2-D gel images using **A)** wide range Immobililine DryStrip pH 3-11NL or **B-E)** four overlapping narrow range IPG strips. Immobililine DryStrip pH 7-11 NL is able to resolve a number of spots that were lost in the basic area of Immobililine DryStrip pH 3-11 NL.

Higher resolution using narrow range Immobililine DryStrip gels

To maximize resolution, four overlapping narrow range IPG strips were used and the spot maps were analyzed using ImageMaster™ 2D Platinum. The analysis was compared to data obtained using the wide range Immobililine DryStrip pH 3-11 NL, and the analyses revealed a 200% increase of detected protein spots. Figure 3 shows the gel images from the narrow and wide range IPG strips. Figure 3E shows that Immobililine DryStrip pH 7-11 NL was able to resolve several spots that were lost in the basic area of Immobililine DryStrip pH 3-11 NL.

Improved IPG Buffer pH 3-10 and pH 3-10 NL increases resolution

IPG Buffer pH 3-10 and pH 3-10 NL have been improved with a new formula. The improved buffer formula was compared with the old formula based on the resulting 2-D spot maps. The quality of the gels was evaluated both by a panel of six experienced 2-D users and by software analysis using ImageMaster 2D Platinum. The panel evaluated several triplicate sets of 2-D spot maps using either improved IPG Buffer or previous IPG Buffer formulations, and based their evaluation on the numbers of spots, intensity, and focusing of the spots. Using ImageMaster 2D Platinum, the 2-D maps were evaluated based on spot intensity, spot matching, and spot numbers. The basic pH area of two 2-D spot maps is shown in Figure 4, while Figure 5 shows the complete 2-D spot maps. Both investigation methods demonstrated that the improved buffers gave additional spots in the high molecular weight area and basic pH area of the gel, as well as spots with increased intensities.

By increasing the intensity of low intensity spots, the probability of detecting low abundant proteins is increased. Note that the improved formula for IPG Buffer pH 3-10 and pH 3-10 NL gives rise to a higher background in the basic pH area of the gel when using silver or Coomassie™ staining. This higher background can be avoided by prolonging the fixation time to 3 hours or up to overnight if needed.

- Sample labeling:* Lysed *E.coli* cells (50 µg) labeled with CyDye™ DIGE Fluor, Cy™5 minimal dye
- Sample buffer:* 8 M urea, 4% CHAPS, 40 mM DTT, and 2% IPG Buffer pH 3-10 NL.
- First dimension IEF:* Immobililine DryStrip gels (pH 3-10 NL; 18 cm) were rehydrated overnight at room temperature with Destreak Rehydration Solution and 0.5% IPG Buffer pH 3-10 NL. The sample was applied by anodic cup loading and focused for 24 kWh.
- Second dimension:* Run on Ettan™ DALTwelve electrophoresis system using 12.5% labcast Laemmli gels.




Fig 4. The basic area of 2-D protein spot maps using Immobililine DryStrip pH 3-10 NL (18 cm) together with **A)** improved IPG Buffer pH 3-10 NL or **B)** the previous recipe for IPG Buffer pH 3-10 NL. The improved IPG Buffer formula results in additional spots, and some of the differences are indicated in the figure.




Fig. 5. 2-D protein spot maps using Immobiline DryStrip pH 3-10 NL (18 cm) together with **A)** improved IPG Buffer pH 3-10 NL or **B)** the previous IPG Buffer pH 3-10 NL formulation. Experimental conditions were identical to the conditions in Figure 4.

Choosing IPG strip and buffer

A comprehensive range of overlapping IPG strips covering narrow (~1 pH unit), medium (3-5 pH units), and wide (~8 pH units) pH ranges are available in five different strip lengths (7, 11, 13, 18, and 24 cm). To simplify gel use and record keeping, each strip is labeled with the pH interval, batch number, individual number, and bar code for use with a reader. Figure 6 shows an overview of the IPG strips available along with the matching IPG Buffers. The combination of IPG strip and IPG Buffer of matching pH interval improves the conductivity distribution across the pH gradient during IEF. The technical specifications for Immobiline DryStrip gels and IPG Buffers are shown in Table 1.

To gain a broad overview of total protein distribution, select wide range IPG strips such as Immobiline DryStrip pH 3-11 NL. Choose shorter strips, up to 13 cm, for fast, cost-effective screening, for a quick overview, or when only the most abundant proteins are of interest (as in pre-fractionated protein samples). The shortest IPG strips

give the fastest results, but their sample load is limited. For more detailed studies, choose medium range IPG strips. To obtain higher resolution, narrow range IPG strips may be used. Use the longer 18 and 24 cm strips for maximal resolution and loading capacity. Longer strips allow detection of more spots and make it easier to select and identify the proteins in the spots, but require longer focusing times. For a more detailed insight into protein distribution, combine medium range Immobiline DryStrip pH 3-7 (or pH 4-7) and pH 6-11 (alternatively, pH 6-9 or pH 7-11). This combination is also useful for preparative applications. Medium range IPG strips in the 24 cm long format deliver both high resolution and high protein loading capacity, both of which improve protein identification and characterization later in the workflow. Deep Purple™, Coomassie, or silver stained spots can be picked from the gel and identified by mass spectrometry.

Table 1. Technical specifications for Immobiline DryStrip gels and IPG Buffer

Immobiline DryStrip	
Gel dimensions	70 × 3 × 0.5 mm, 110 × 3 × 0.5 mm, 130 × 3 × 0.5 mm, 180 × 3 × 0.5 mm, 240 × 3 × 0.5 mm
Gel matrix	Polyacrylamide T = 4%, C = 3%
Gel backing	Polyester film
Storage	-20°C
Shelf life	18 months
IPG Buffer	
Content	Specialized carrier ampholytes in aqueous solution
Volume	1 ml
Storage	4°C to 8°C
Shelf life	36 months

	IPG Strips					IPG Buffers							
	Strip length					pH range							
	24 cm	18 cm	13 cm	11 cm	7 cm	3.5-5.0	5.5-6.7	4-7	6-11	7-11 NL	3-10 NL	3-10	3-11 NL
Narrow													
3.5-4.5	×					•							
5.3-6.5	×	×	×	×	×		•						
6.2-7.5	×	×	×	×	×				•				
Medium													
3-5.6 NL	×	×	×	×	×	•							
3-7 NL	×							•					
4-7	×	×	×	×	×			•					
6-9	×	×							•				
6-11		×	×	×	×				•				
7-11 NL	×	×	×	×	×					•			
Wide													
3-10	×	×	×	×	×							•	
3-11 NL	×	×	×	×	×								•
3-10 NL	×	×	×	×	×						•		

Fig 6. A comprehensive range of overlapping IPG strips covering narrow, medium, and wide pH ranges are available in several different strip lengths. Using IPG strips together with matching IPG Buffer improves the conductivity distribution across the pH gradient during IEF.

IPGbox

Immobiline DryStrip gel packs are shipped dry for stability. IPGbox provides a convenient method for rehydrating up to twelve precast IPG strips (7 to 24 cm) at a time. Individual slots in the Reswell Trays allow rehydration of individual IPG strips in a minimum volume of solution, and no oil is required. IPGbox can accommodate all IPG strip sizes including 7, 11, 13, 18, and 24 cm (Fig 7). IPGbox is supplied with a complete set of accessories to rehydrate 10 full rounds of 12 strips each. To avoid redox-related streaking, IPG strips should be rehydrated with DeStreak Rehydration Solution.



Fig 7. IPGbox is specially designed for rehydrating Immobiline DryStrip gels in a minimum volume of solution.

Synchronized solutions for better 2-D electrophoresis - improved 2-D DIGE

2-D electrophoresis is a powerful technology for comparing complex protein mixtures from biological samples in proteomics research. With the aim of simplifying the process and improving quantitation, GE Healthcare has made significant improvements across the entire workflow, from sample prep and separation to data analysis. 2-D DIGE

(Difference Gel Electrophoresis) allows you to analyze differences in low abundance proteins with unparalleled statistical confidence. Many of our 2-D electrophoresis products improve data quality from traditional 2-D experiments as well, and the products can eventually be integrated into a complete solution to maximally improve your 2-D results.



Sample Preparation	Protein Labeling	Protein Separation	Image Acquisition	Image Analysis	Validation
2-D Protein Extraction Buffer	DIGE Trial Pack CyDye™ DIGE Fluor 2 nmol kit CyDye DIGE 400 nmol sat dye CyDye DIGE repackaging	IPGbox Immobiline™ DryStrip gels IPG Buffer DIGE Gel and DIGE Buffer Kit	Ettan™ DIGE Imager Typhoon™ Variable Mode Imager	DeCyder™ 2-D software v7.0	ECL Plex™ ImageQuant™ TL 7.0 with ImageQuant TL SecurITy

For more information, visit www.gelifesciences.com/2DE

Ordering information

Product	Quantity	Code No	Related products	Quantity	Code No
IPG Buffer pH 3.5-5.0	1 ml	17-6002-02	IPGbox	1 IPGbox + 1 IPGbox Kit	28-9334-65
IPG Buffer pH 5.5-6.7	1 ml	17-6002-06			
IPG Buffer pH 4-7	1 ml	17-6000-86	IPGbox Kit	10 Reswell Trays + 1 IPGbox Insert	28-9334-92
IPG Buffer pH 6-11	1 ml	17-6001-78			
IPG Buffer pH 7-11 NL	1 ml	17-6004-39			
IPG Buffer pH 3-10	1 ml	17-6000-87	DeStreak Rehydration Solution	5 × 3 ml	17-6003-19
IPG Buffer pH 3-10 NL	1 ml	17-6000-88	DeStreak Reagent	1 ml	17-6003-18
IPG Buffer pH 3-11 NL	1 ml	17-6004-40	Ettan IPGphor™ 3 Isoelectric Focusing Unit	1	11-0033-64
Related literature		Code No	Ettan IPGphor Manifold, Complete	1	80-6498-38
2-D Electrophoresis: Principles and Methods, Handbook		80-6429-60	Ettan IPGphor Manifold, Light Complete	1	11-0026-88
Ettan DIGE System User Manual		18-1173-17			

Immobiline DryStrip gels¹

	Code No					
Narrow range	pH 3.5-4.5	pH 5.3-6.5	pH 6.2-7.5			
24 cm	17-6002-38	17-6003-62	17-6003-67			
18 cm	-	17-6003-61	17-6003-66			
13 cm	-	17-6003-60	17-6003-65			
11 cm	-	17-6003-59	17-6003-64			
7 cm	-	17-6003-58	17-6003-63			
Medium range	pH 3-5.6 NL	pH 3-7 NL	pH 4-7	pH 6-9	pH 6-11	pH 7-11 NL
24 cm	17-6003-57	17-6002-43	17-6002-46	17-6002-47	-	17-6003-72
18 cm	17-6003-56	-	17-1233-01	17-6001-88	17-6001-97	17-6003-71
13 cm	17-6003-55	-	17-6001-13	-	17-6001-96	17-6003-70
11 cm	17-6003-54	-	18-1016-60	-	17-6001-95	17-6003-69
7 cm	17-6003-53	-	17-6001-10	-	17-6001-94	17-6003-68
Wide range	pH 3-10	pH 3-10 NL	pH 3-11 NL			
24 cm	17-6002-44	17-6002-45	17-6003-77			
18 cm	17-1234-01	17-1235-01	17-6003-76			
13 cm	17-6001-14	17-6001-15	17-6003-75			
11 cm	18-1016-61	-	17-6003-74			
7 cm	17-6001-11	17-6001-12	17-6003-73			

¹ Pack size is 12 strips/pack

For local office contact information, visit
www.gelifesciences.com/contact

www.gelifesciences.com/2DE

GE Healthcare Bio-Sciences AB
Björkgatan 30
751 84 Uppsala
Sweden



GE, imagination at work, and GE monogram are trademarks of General Electric Company.

Cy, CyDye, DeCyder, Deep Purple, ECL Plex, Ettan, ImageMaster, ImageQuant, Immobiline, IPGphor, and Typhoon are trademarks of GE Healthcare companies.

2-D Fluorescence Difference Gel Electrophoresis (2-D DIGE) technology is covered by US patent numbers 6,043,025, 6,127,134 and 6,426,190 and equivalent patents and patent applications in other countries and exclusively licensed from Carnegie Mellon University. CyDye: this product or portions thereof is manufactured under an exclusive license from Carnegie Mellon University under US patent numbers 5,569,587, 5,627,027 and equivalent patents in other countries. The purchase of CyDye DIGE Fluors includes a limited license to use the CyDye DIGE Fluors for internal research and development, but not for any commercial purposes. A license to use the CyDye DIGE Fluors for commercial purposes is subject to a separate license agreement with GE Healthcare.

CyDye: This product or portions thereof is manufactured under an exclusive license from Carnegie Mellon University under US patent number 5,268,486 and equivalent patents in the US and other countries.

The purchase of CyDye products includes a limited license to use the CyDye products for internal research and development but not for any commercial purposes. A license to use the CyDye products for commercial purposes is subject to a separate license agreement with GE Healthcare. Commercial use shall include:

1. Sale, lease, license or other transfer of the material or any material derived or produced from it.
2. Sale, lease, license or other grant of rights to use this material or any material derived or produced from it.
3. Use of this material to perform services for a fee for third parties, including contract research and drug screening.

If you require a commercial license to use this material and do not have one, return this material unopened to GE Healthcare Bio-Sciences AB, Björkgatan 30, SE-751 84 Uppsala, Sweden and any money paid for the material will be refunded.

DeCyder: This release of DeCyder (software) is provided by GE Healthcare to the customer under a nonexclusive license and is subject to terms and conditions set out in the 2-D Differential Gel Electrophoresis Technology Access Agreement. Customer has no rights to copy or duplicate or amend the Software without the prior written approval of GE Healthcare.

Deep Purple Total Protein Stain is exclusively licensed to GE Healthcare from Fluorotronics Pty Ltd. Deep Purple Total Protein Stain may only be used for applications in life science research. Deep Purple is covered under a granted patent in New Zealand entitled "Fluorescent Compounds", patent number 522291 and equivalent patents and patent applications in other countries.

DIGE Gel and DIGE Buffer Kit: The buffer system in this gel and buffer kit is covered by patent application WO9616724 granted in US, EP and JP.

This version of ImageMaster has been developed by the Swiss Institute of Bioinformatics in collaboration with GeneBio and GE Healthcare.

All third party trademarks are the property of their respective owners.

© 2009 General Electric Company—All rights reserved.

First published Mar. 2009

All goods and services are sold subject to the terms and conditions of sale of the company within GE Healthcare which supplies them. A copy of these terms and conditions is available on request. Contact your local GE Healthcare representative for the most current information.

GE Healthcare UK Limited Amersham Place
Little Chalfont
Buckinghamshire, HP7 9NA
UK

GE Healthcare Europe, GmbH
Munzinger Strasse 5
D-79111 Freiburg
Germany

GE Healthcare Bio-Sciences Corp.
800 Centennial Avenue, P.O. Box 1327
Piscataway, NJ 08855-1327
USA

GE Healthcare Bio-Sciences KK
Sanken Bldg., 3-25-1, Hyakunincho
Shinjuku-ku, Tokyo 169-0073
Japan