

GE Healthcare



Combine 4 instruments into 1 without compromising sensitivity and dynamic range

Amersham™ Typhoon™ Biomolecular Imager



Detect low levels of proteins

The ability to distinguish subtle differences in expression among low amounts of biomolecules in complex mixtures adds power to many areas of biological research. The highly versatile Amersham Typhoon laser scanner enables you to generate data of the highest quality through linearity of signal response, quantitative accuracy, and extremely low limits of detection. The completely customizable Typhoon supports phosphorimaging, 2-D DIGE imaging, red/green/blue (RGB) and near IR fluorescence as well as sensitive and accurate quantitation of proteins.

Enhance your sensitivity

- Improved PMTs with 5 Laser Configuration Option
- Multiplex detection of proteins with minimal crosstalk




Fig. 1. Phosphorylase B was labeled with CyDye™ DIGE fluor Cy5 minimal dye and separated using the Amersham WB electrophoresis gel. The gel was imaged with Amersham Typhoon using normal scan speed. The detection limit was 5.6 pg and the linear dynamic range (DR) was 4.8 orders of magnitude.

Enhance your dynamic range

- Broad linear dynamic range greater than 5 orders of magnitude




Fig. 2. Phosphorylase B was labeled with CyDye DIGE fluor Cy3 minimal dye and separated using the Amersham WB electrophoresis gel. The gel was imaged with Amersham Typhoon using normal scan speed. The detection limit was 22 pg and the linear dynamic range (DR) was 4.2 orders of magnitude.



Three preconfigured models available:

	Phosphor Imaging	Densitometry (OD)	RGB Fluorescence	IR1 & IR2 Fluorescence
Amersham Typhoon IP	✓	○	○	○
Amersham Typhoon RGB	✓	✓	✓	○
Amersham Typhoon 5	✓	✓	✓	✓

4 Instruments into 1

The Amersham Typhoon is an upgradeable laser scanner, completely customizable for your specific research needs. The system combines 4 instruments into 1, which allows for several imaging modes at any time, including near infrared (NIR) fluorescence, RGB fluorescence, Phosphor imaging, and densitometry of colorimetrically-stained samples (e.g., Coomassie™ blue and silver stain).

Enhance your flexibility

- Select the lasers you need for blue, green, red, or NIR detection
- Phosphorimaging for radiolabeled target detection and enhanced sensitivity




Fig. 3. Antibody conjugated with IRDye™ 680 was separated using the Amersham WB electrophoresis gel. To reduce noise, the gel was imaged with Amersham Typhoon using slow scan speed. The detection limit was 3 pg and the linear dynamic range (DR) was 5.2 orders of magnitude.




Fig. 4. Antibody conjugated with IRDye™ 800 was separated using the Amersham WB electrophoresis gel. To reduce noise, the gel was imaged with Amersham Typhoon using slow scan speed. The detection limit was 3 pg and the linear dynamic range (DR) was 5.2 orders of magnitude.




Fig. 5. Autoradiography images of rat injected with ^{14}C glucose. The magnified area shows part of the spine.




Fig. 6. Scanned image of a ^{14}C autoradiographic standard using Amersham Typhoon. The linear dynamic range (DR) was 5.3 orders of magnitude.

Amersham Eraser

Amersham Eraser is a dedicated instrument for erasing imaging plates (storage phosphor screens), saving time and cost compared to using film. The LED light source provides an environmentally friendly option, with less energy consumption and longer life time than conventional light sources, while the drawer design with touch controls provides an easy to use interface.



Ordering information

System	Quantity	Product code
Amersham Typhoon 5	1	29187191
Amersham Typhoon RGB	1	29187193
Amersham Typhoon IP	1	29187194
Amersham Eraser	1	29187190
Accessory Cabinet AmTyphoon	1	29191637

GE Healthcare Bio-Sciences Corp.
100 Results Way
Marlborough
MA 01752
USA
www.gelifesciences.com/typhoon

www.gelifesciences.com

GE, GE monogram, Amersham, Cy, CyDye, and Typhoon are trademarks of General Electric Company.

IRDye is a trademark of LI-Cor Biosciences, Inc.

Coomassie is a trademark of Thermo Fisher Scientific, Inc.

All other third party trademarks are the property of their respective owner.

© 2016 General Electric Company.

First published Sept. 2016

All goods and services are sold subject to the terms and conditions of sale of the company within GE Healthcare which supplies them. A copy of these terms and conditions is available on request. Contact your local GE Healthcare representative for the most current information.

GE Healthcare Bio-Sciences AB, Björkgatan 30, SE-751 84 Uppsala, Sweden

GE Healthcare UK Limited, Amersham Place, Little Chalfont, Buckinghamshire, HP7 9NA, UK

GE Healthcare Dharmacon, Inc., 2650 Crescent Dr., Lafayette, CO 80026, USA

HyClone Laboratories, Inc., 925 W 1800 S, Logan, UT 84321, USA

GE Healthcare Europe GmbH, Munzinger Strasse 5, D-79111 Freiburg, Germany

GE Healthcare Japan Corporation, Sanken Bldg. 3-25-1, Hyakunincho, Shinjuku-ku, Tokyo 169-0073, Japan

For local office contact information, visit: www.gelifesciences.com/contact

29222858AA 09/2016

